



The coil connector socket allows a treatment coil to safely and securely be locked into place.



Magnus Compact Brief User Guide

The MAINS power rocker switch turns the unit on "I" and off "O". MAINS socket to plug the medical grade MAINS power cord into. Always plug the power cord into a surge protected MAINS outlet.

The pulse rate/intensity dial increases the intensity of the electromagnetic pulse as it decreases the pulse rate. The pulse intensity can be adjusted from low (1) to high (10).

Preset 5-minute and 10-minute pushbuttons begin treatment cycles. The blue light turns on when the cycle begins

DIRECTIONS FOR USE

Getting Started

- ⚠ Always connect a surge protector between the unit and the MAINS.
- ⚠ Connect the power cord into the MAINS socket on the inside of the unit.
- ⚠ Always lock a treatment coil connector into the connector socket on the panel of the unit before pressing the MAINS power switch on "I". The unit will not start if a treatment coil is not properly connected.

Inserting/Changing a Treatment Coil

Do not unplug or plug in a treatment coil during the treatment cycle.

Insert the treatment coil connector by lining up the flat top of the connector with the coil socket then push it in until you hear a click. The treatment coil is then locked into place.

Remove the treatment coil connector by depressing the connector spring and pull the coil connector straight out of the socket.

Beginning Treatment

Press the MAINS power rocker switch to turn the unit on "I".

Press the power switch again to turn the unit off "O".

Place the treatment coil on or around the desired treatment area. The closer the coil is to the treatment area, the more effective the treatment will be. Users can remain fully clothed and no direct contact between the coil and the skin is necessary for the treatment to be effective.

Make sure the pulse rate/intensity dial is on 1 before beginning treatment on any area.

Press the 5 or 10-minute start pushbutton to begin treatment. A blue light turns on when the pushbutton is pressed. The blue light turns off and the unit beeps when the treatment ends. Press one of the start pushbuttons again to begin another treatment cycle.

After beginning treatment, raise the intensity to a comfortable level for treatment by turning the dial towards 10. Lower the intensity back to 1 when moving the coil to a different area.

Pausing/Resetting/Ending the Treatment

Pressing the illuminated start pushbutton after the treatment cycle has begun places the unit on pause and the blue light flashes.

Pressing the start pushbutton once more resumes the treatment cycle.

Holding the start pushbutton down until it beeps, stops the treatment cycle and the blue light turns off.

A single beep sound signifies the end of the treatment cycle.

CONTRAINDICATIONS

DO NOT USE IF:

- You are, or may be pregnant,
- You are receiving chemotherapy,
- You have cancer or malignant tumors,
- You have had surgery in the past 24 hours,
- You have hemorrhagic tendencies, Purpura or Hemophilia,
- You have major metabolic diseases uncontrolled by medication (HIV, ulcers...).

DO NOT USE within 1 ft. of pacemakers, defibrillators or any other implanted electronic devices.

DO NOT USE within 1 ft. of metallic implants (stents, pins, rods, or screws).

DO NOT USE soon after taking any medication.

WARNINGS & PRECAUTIONS

-  Attention: Dangerous Voltage inside, DO NOT tamper or open the control unit.
-  Always connect a surge protector between the unit and the MAINS.
-  Do not use the unit near credit cards, security access cards, car keys, hearing aids, watches, jewelry and metal chains, cell phones, IPODS, laptops, remote controls, or any other electronic media. The electromagnetic fields may disrupt their functioning and/or demagnetize them.
-  Do not use in wet environments. Do not expose to moisture or liquids of any kind.
-  Do not block fan vents.
-  Do not apply stimulation while operating any machinery or during any activity in which involuntary muscle contractions may put the user at undue risk of injury.
-  Keep the device out of the reach of children.
-  Caution: Users with low blood pressure may feel temporarily dizzy when first standing up.
-  Caution: Users with inactive lifestyles should not use the unit for more than one 5-minute treatment cycle per area for the first few times to prevent sore muscles.
-  Caution: Users with heart problems should consult their physician before using.
-  Caution: Users with suspected or diagnosed epilepsy should consult their physician before using.
-  Caution: Do not place the treatment coil over a suture line within 3 days after surgery.
-  Caution: Users taking pain, anxiety, depression or any other medication should be carefully monitored when using the unit as medication effectiveness may be intensified.
-  Caution: Applying directly over the menstruating uterus may cause increased bleeding.
-  Caution: Do not apply at high intensity over areas of the skin or body that lack normal sensation.
-  Caution: Metal objects will heat with prolonged exposure to electromagnetic fields. Use the one 5-minute treatment cycle on low power followed by a five minute pause within 1 ft. of any implanted metal objects.

ADVERSE REACTIONS

There are no known or reported adverse reactions with the use of the unit or any other electromagnetic stimulation devices.

Detoxification may occur after treatment. Drink plenty of water after treatment.

Overuse may cause muscle soreness.

Users with low blood pressure may experience temporary dizziness when first standing up.

In case of adverse reaction, stop use and consult a physician.